

# Pancreatic Stents

Single Use

## Pancreatic Stents

- Well Balanced Wall Diameter  
*Maximum Inner Lumen with Highest Form Stability & Kink Resistance*
- Tapered & Round Ends  
*Atraumatic Implantation*
- Multiple Holes  
*Perfect Drainage of Obstructed Pancreatic Ducts*
- Radiopaque Material  
*Exact Control of Placement*




## Pancreatic Stents

Medi-Globe offers the pancreatic stents as single component (catalog no. without "S") and additionally as integrated part of a stent set.

A stent set comprises of:


- Stent
- Pusher



### Pancreatic Stent Set

Straight with Radial Holes, Distal and Duodenal Flaps


Catalog No.	Stent		Pusher	
	Size / French	Length	Size / French	Length
GBS-10-05-003S	5.0	3 cm	5.0	220 cm
GBS-10-05-005S	5.0	5 cm	5.0	220 cm
GBS-10-05-007S	5.0	7 cm	5.0	220 cm
GBS-10-05-009S	5.0	9 cm	5.0	220 cm
GBS-10-05-011S	5.0	11 cm	5.0	220 cm
GBS-10-07-003S	7.0	3 cm	7.0	220 cm
GBS-10-07-005S	7.0	5 cm	7.0	220 cm
GBS-10-07-007S	7.0	7 cm	7.0	220 cm
GBS-10-07-008S	7.0	8 cm	7.0	220 cm
GBS-10-07-009S	7.0	9 cm	7.0	220 cm
GBS-10-07-011S	7.0	11 cm	7.0	220 cm
GBS-10-07-016S	7.0	16 cm	7.0	220 cm
GBS-10-08-005S	8.5	5 cm	8.5	220 cm
GBS-10-08-007S	8.5	7 cm	8.5	220 cm
GBS-10-08-008S	8.5	8 cm	8.5	220 cm
GBS-10-08-009S	8.5	9 cm	8.5	220 cm
GBS-10-08-011S	8.5	11 cm	8.5	220 cm
GBS-10-08-016S	8.5	16 cm	8.5	220 cm



### Pancreatic Stent Set

Straight with Radial Holes, One Duodenal Flap


Catalog No.	Stent		Pusher	
	Size / French	Length	Size / French	Length
GBS-20-05-003S	5.0	3 cm	5.0	220 cm
GBS-20-05-005S	5.0	5 cm	5.0	220 cm
GBS-20-05-007S	5.0	7 cm	5.0	220 cm



### Pancreatic Stent Set

Single Pigtail with Multiple-Flaps, Multiple Drainage Holes, Position Marker

Catalog No.	Stent		Pusher	
	Size / French	Length	Size / French	Length
GBS-08-05-003S	5.0	3 cm	5.0	220 cm
GBS-08-05-004S	5.0	4 cm	5.0	220 cm
GBS-08-05-005S	5.0	5 cm	5.0	220 cm
GBS-08-05-006S	5.0	6 cm	5.0	220 cm
GBS-08-05-007S	5.0	7 cm	5.0	220 cm
GBS-08-05-008S	5.0	8 cm	5.0	220 cm
GBS-08-05-009S	5.0	9 cm	5.0	220 cm
GBS-08-05-010S	5.0	10 cm	5.0	220 cm
GBS-08-05-012S	5.0	12 cm	5.0	220 cm
GBS-08-07-002S	7.0	2 cm	7.0	220 cm
GBS-08-07-004S	7.0	4 cm	7.0	220 cm
GBS-08-07-006S	7.0	6 cm	7.0	220 cm
GBS-08-07-009S	7.0	9 cm	7.0	220 cm
GBS-08-07-012S	7.0	12 cm	7.0	220 cm
GBS-08-08-006S	8.5	6 cm	8.5	220 cm
GBS-08-08-009S	8.5	9 cm	8.5	220 cm
GBS-08-08-012S	8.5	12 cm	8.5	220 cm
GBS-08-10-006S	10.0	6 cm	10.0	220 cm
GBS-08-10-009S	10.0	9 cm	10.0	220 cm
GBS-08-10-012S	10.0	12 cm	10.0	220 cm
GBS-08-11-006S	11.5	6 cm	11.5	220 cm
GBS-08-11-009S	11.5	9 cm	11.5	220 cm
GBS-08-11-012S	11.5	12 cm	11.5	220 cm




### Pancreatic Stent Set

Single Pigtail without Flap, Multiple Drainage Holes, Position Marker

Catalog No.	Stent		Pusher	
	Size / French	Length	Size / French	Length
GBS-18-05-003S	5.0	3 cm	5.0	220 cm
GBS-18-05-004S	5.0	4 cm	5.0	220 cm
GBS-18-05-005S	5.0	5 cm	5.0	220 cm
GBS-18-05-006S	5.0	6 cm	5.0	220 cm
GBS-18-05-007S	5.0	7 cm	5.0	220 cm

# Pancreatic Stents

Single Use



## Pancreatic Stent Set S - Shape

Catalog No.	Stent		Pusher	
	French Size	Length	French Size	Length
GBS-09-05-005S	5.0	5 cm	5.0	220 cm
GBS-09-05-007S	5.0	7 cm	5.0	220 cm
GBS-09-05-009S	5.0	9 cm	5.0	220 cm
GBS-09-05-012S	5.0	12 cm	5.0	220 cm
GBS-09-05-015S	5.0	15 cm	5.0	220 cm
GBS-09-07-005S	7.0	5 cm	7.0	220 cm
GBS-09-07-007S	7.0	7 cm	7.0	220 cm
GBS-09-07-009S	7.0	9 cm	7.0	220 cm
GBS-09-07-012S	7.0	12 cm	7.0	220 cm
GBS-09-07-015S	7.0	15 cm	7.0	220 cm

## Pushers

Pushers follow for each diameter a defined color concept to differentiate the related French sizes.



## Pusher with 1 Metal Ring

Catalog No.	French Size	Color Concept	Length
GBS-42-05-220	5.0	grey	220 cm
GBS-42-07-160	7.0	yellow	160 cm
GBS-42-07-220	7.0	yellow	220 cm
GBS-42-08-160	8.5	orange	160 cm
GBS-42-08-180	8.5	orange	180 cm
GBS-42-08-220	8.5	orange	220 cm
GBS-42-10-160	10.0	red	160 cm
GBS-42-10-180	10.0	red	180 cm
GBS-42-10-220	10.0	red	220 cm
GBS-42-11-160	11.5	blue	160 cm
GBS-42-11-180	11.5	blue	180 cm
GBS-42-11-220	11.5	blue	220 cm

Distributed by: