



aerstent®

trachea bronchus – wire woven design

Trachea – Bronchus stents

Leufen Medical offers a huge variety of self-expandable tracheal and bronchial stents with diameters from 10 to 20 mm and lengths from 20 to 80 mm. The spherical end design and the strong radial force keep the stents in position and avoid migration.

NEW To reduce the risk of mucostasis and the colonization of pathogenic germs the stent is hydrophilic internal coated with **Lubriteq®**.

Practical application

All tracheal and bronchial stents are preloaded on the **new** reload introducer system. For correction of the stent position the stent can be reloaded into the introducer system during the procedure. The excellent radiopacity makes the stents clearly visible to x-ray.

Made in Germany

With our new product line **aerstent®** you can expect a high quality product made in Germany which meets all of your ambitious requirements fully.

There are three sectors in which Leufen Medical is acting as a role model and stent pioneer

- innovative development
- precise manufacture
- patient specific product distribution

Order information you find on the reverse side.
Revision status September 2010

aerstent[®]

trachea bronchus – wire woven design

- **NEW** stent with complete covering
- **NEW** hydrophilic internal coated **Lubriteq[®]**
- Atraumatic, rounded ends
- Memory effect of nitinol material – high radial force and flexibility
- **NEW** reload introducer
 - Length 600 mm, OTW
 - Ø 6 mm / 18 French (stents from 10 to 14 mm)
 - Ø 8 mm / 24 French (stents from 16 to 20 mm)
 - Accepts guide wire 0.035 inch

trachea bronchus – complete covered

REF	Length	Length Cover	Ø Middle	Ø Ends
503-10-030	30	30	10	12
503-12-020	20	20	12	14
503-12-040	40	40	12	14
503-14-020	20	20	14	16
503-14-030	30	30	14	16
503-14-040	40	40	14	16
503-16-030	30	30	16	18
503-16-040	40	40	16	18
503-16-050	50	50	16	18
503-16-060	60	60	16	18
503-16-080	80	80	16	18
503-18-030	30	30	18	20
503-18-050	50	50	18	20
503-18-060	60	60	18	20
503-18-080	80	80	18	20
503-20-040	40	40	20	22
503-20-060	60	60	20	22
503-20-080	80	80	20	22
503-22-060	60	60	22	24

Leufen Medical GmbH
 Kapellenstr. 82
 D - 52066 Aachen / Germany
 T +49 (0)2 41 - 997 10 10
 F +49 (0)2 41 - 997 10 11
 contact@leufen-medical.de

