



The anastomosis colon and rectum stents are preloaded on the OTW (over-the-wire) delivery systems. Their specific end designs have the perfect fit for the anatomic and peristaltic conditions in the colon and the rectum.

The anastomosis colon and rectum stents are preloaded on the reload delivery system. For correction of the stent position the stent can be reloaded into the delivery system during the procedure. The excellent radiopacity tantalum marker makes the stents clearly visible in the fluoroscopy.

- Stent with complete silicone covering
- Atraumatic ends
- High tensile extraction threads on both ends
- Special tulip design prevent migration
- Memory effect of nitinol material – high radial force
- Radiopaque tantalum markers
- Reload delivery system

Length 1100 mm, OTW  
 Ø 8 mm / 24 French  
 Accepts guide wire 0.035 inch

aixstent · colon rectum anastomosis · complete silicone covered

REF	Length	Length Cover	Ø Middle	Ø Ends
403-20-060	60	60	20	26
403-30-060	60	60	30	36