

Cytology Brushes

Single Use

Cytology Brushes

- Brush with Metal Tip
Efficient Tube Closure to Secure Smear Yield
- Set Versions including Brush and Scissor
Safe and User-friendly Specimen Handling
- Sophexa Easy Brush
Guide Wire Compatible Cytology Brush for Smear Acquisition from the Bile Duct



Cytology Brushes

Single Use

Bronchoscopy

Cytology Brush – with Conical Metal Tip and Ergonomic Handle

| Catalog No. | Description | Working Length | Sheath Outer Diameter | Minimum Accessory Channel Diameter |
|---------------|------------------------------------|----------------|-----------------------|------------------------------------|
| GCB-01-18-120 | Cytology Brush | 120 cm | 1.8 mm | 2.0 mm |
| GCB-01-24-120 | Cytology Brush | 120 cm | 2.45 mm | 2.8 mm |
| GCB-02-18-120 | Set incl. Cytology Brush & Scissor | 120 cm | 1.8 mm | 2.0 mm |
| GCB-02-24-120 | Set incl. Cytology Brush & Scissor | 120 cm | 2.45 mm | 2.8 mm |

Gastroscopy

Cytology Brush – with Conical Metal Tip and Ergonomic Handle

| Catalog No. | Description | Working Length | Sheath Outer Diameter | Minimum Accessory Channel Diameter |
|---------------|------------------------------------|----------------|-----------------------|------------------------------------|
| GCB-01-18-180 | Cytology Brush | 180 cm | 1.8 mm | 2.0 mm |
| GCB-01-24-180 | Cytology Brush | 180 cm | 2.45 mm | 2.8 mm |
| GCB-02-18-180 | Set incl. Cytology Brush & Scissor | 180 cm | 1.8 mm | 2.0 mm |
| GCB-02-24-180 | Set incl. Cytology Brush & Scissor | 180 cm | 2.45 mm | 2.8 mm |

Colonoscopy

Cytology Brush – with Conical Metal Tip and Ergonomic Handle

| Catalog No. | Description | Working Length | Sheath Outer Diameter | Minimum Accessory Channel Diameter |
|---------------|------------------------------------|----------------|-----------------------|------------------------------------|
| GCB-01-18-230 | Cytology Brush | 230 cm | 1.8 mm | 2.0 mm |
| GCB-01-24-230 | Cytology Brush | 230 cm | 2.45 mm | 2.8 mm |
| GCB-02-18-230 | Set incl. Cytology Brush & Scissor | 230 cm | 1.8 mm | 2.0 mm |
| GCB-02-24-230 | Set incl. Cytology Brush & Scissor | 230 cm | 2.45 mm | 2.8 mm |

Sophexa Easy Brush

Guide Wire Compatible Cytology Brush for Smear Acquisition mainly from the Bile Duct.

Cytology Brush – with Ellipsoid Tip and Ergonomic Handle

| Catalog No. | Description | Working Length | Sheath Outer Diameter | Minimum Accessory Channel Diameter |
|---------------|----------------|----------------|-----------------------|------------------------------------|
| GCB-03-24-230 | Cytology Brush | 230 cm | 2.45 mm | 2.8 mm |