

# Nasal Biliary Drainages

Single Use

## Nasal Biliary Drainages

- Various Types and Sizes Available  
*Fit for All Anatomic Cases*
- Well Balanced Wall Diameter  
*Maximum Inner Lumen with Highest Form Stability & Kink Resistance*
- Tapered & Round Ends  
*Atraumatic Implantation*



# Nasal Biliary Drainages

Single Use

## Nasal Biliary Drainages Sets

### Nasal Biliary Drainages Sets

Set Components: Drainage Catheter, Nasal Tube, Drainage Connecting Tube, Injection Connector and Teflon Coated Stainless Steel Guide Wire

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube	Drainage Connecting Tube	Guide Wire	
		French Size	Length	Length	Length	Diameter	Length
GND-01-05-290S	Pigtail	5.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-01-07-290S	Pigtail	7.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-01-10-290S	Pigtail	10.0	290 cm	50 cm	135 cm	.052" (1.32 mm)	400 cm
GND-02-05-290S	Standard	5.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-02-07-290S	Standard	7.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-02-10-290S	Standard	10.0	290 cm	50 cm	135 cm	.052" (1.32 mm)	400 cm
GND-03-05-290S	Standard with Pigtail	5.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-03-07-290S	Standard with Pigtail	7.0	290 cm	50 cm	135 cm	.035" (0.89 mm)	400 cm
GND-03-10-290S	Standard with Pigtail	10.0	290 cm	50 cm	135 cm	.052" (1.32 mm)	400 cm


All Nasal Biliary Drainage Catheters are also available individually (Catalog set no. without "S").

## Individual Components


### Nasal Tube

Catalog No.	French Size	Length	Product
GND-41-05-050	5.0	50 cm	
GND-41-07-050	7.0	50 cm	
GND-41-10-050	10.0	50 cm	

### Injection Connector

Catalog No.	French Size	Product
GND-43-05-000	5.0	
GND-43-07-000	7.0	
GND-43-10-000	10.0	

### Drainage Connecting Tube

Catalog No.	Length	Product
GND-42-03-135	135 cm	

Ver. 3.0 (2010-09)