

WeightRight® Mercury-Free Bougies



Heavy-weight control and convenience combined.

Our tungsten-filled bougie mimics the weight and feel of mercury-filled bougies without the concerns of mercury usage and disposal. It's the weight you expect, and the convenience you demand.

Benefits

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Mercury-free and latex-free

Features

- 23 sizes: 16 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini set or individual pieces
- Three-year warranty on product and workmanship
- Easy disposal
- No charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)




Store and transport WeightRight® bougies in our storage case.

WeightRight® Dilatation Systems


Set	Part No.	Contents
Maloney (tapered) 23-piece full set	1208-00	23 bougies, 16-60 FR
Maloney (tapered) 10-piece mini set	1208-10	10 bougies, 36-54 FR
Hurst (blunt) 23-piece full set	1209-00	23 bougies, 16-60 FR
Hurst (blunt) 10-piece mini set	1209-10	10 bougies, 36-54 FR
Bougie Storage Case	1210-01	Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H



Our storage case allows easy identification of bougie size. Removing plexiglass window provides easy access.



Hang bougie vertically using the hole in the end cap.



Top: Maloney tip. Bottom: Hurst tip.

Hurst Bougie

Maloney Bougie

American Depth Markings

European Depth Markings



WeightRight® Mercury-free Bougies

FR (French) Size	(mm)	Maloney (tapered) Tip Part No.	Hurst (blunt) Tip Part No.
16	5.3	1208-16	1209-16
18	6	1208-18	1209-18
20	6.6	1208-20	1209-20
22	7.4	1208-22	1209-22
24	8	1208-24	1209-24
26	8.6	1208-26	1209-26
28	9.3	1208-28	1209-28
30	10	1208-30	1209-30
32	10.7	1208-32	1209-32
34	11.2	1208-34	1209-34
36	12	1208-36	1209-36
38	12.5	1208-38	1209-38
40	13	1208-40	1209-40
42	14	1208-42	1209-42
44	14.5	1208-44	1209-44
46	15.2	1208-46	1209-46
48	16	1208-48	1209-48
50	16.5	1208-50	1209-50
52	17.2	1208-52	1209-52
54	18	1208-54	1209-54
56	18.5	1208-56	1209-56
58	19.2	1208-58	1209-58
60	20	1208-60	1209-60