




optiLube™

Lubricating Jelly

water soluble for
general lubrication



www.optimummedical.co.uk



optimum medical solutions



optiLube™

Lubricating Jelly

water soluble for general lubrication

OptiLube™ is a Class IIa sterile lubricating jelly available in a comprehensive range of formats; tubes, syringes and sachets. The lubricant jelly has been developed by clinicians for clinicians. OptiLube™ is designed to meet the specific needs of a range of different procedures, making it the gel of choice, in over 80% of UK hospitals and 36 countries.

OptiLube™:

- Water soluble, pH balanced universal sterile lubricating jelly
- Non-sticky, non-greasy, transparent and petroleum free
- Will not damage rubber or metal and is safe to use on human tissue
- Adheres well to gloves/instruments

OptiLube™ is an essential universal jelly and can be used for a large range of clinical procedures. Some of the possible uses for OptiLube™ include:

Procedure	Insertion/use of	Hospital Department
PAP tests	Speculums	Gynaecology
Surgical procedures	Laryngeal Mask Airways	Operating theatre
Catherisation	Catheters	Urology
Endoscopy	Endoscopes	Endoscopy
Sigmoidoscopy	Sigmoidoscopes	Endoscopy
Cystoscopy	Cystoscopes	Urology
Enemas/douches	Enemas/douches	Proctology
Sonographys (pregnancy, breast imaging etc)	Sonographic Machines	Obstetrics
X-rays (Proctograms)	Photoelectron Spectroscopy	Radiology

Sachets

2.7g e 5g e

OMS Code	Product Description	Brand	Inner Pack	Outer Pack
1123	2.7g Sachet	OptiLube™	Box of 144	Case of 1440
1120	5g Sachet	OptiLube™	Box of 150	Case of 1500

OptiLube™ sachets promote advanced infection control procedures; the individual use sachets eliminate the risk of patient to patient contamination. Available in 2.7g or 5g sachets.

- Individual use ensuring infection control, one sachet per patient
- Easy Tear



Syringes

6ml e 11ml e

OMS Code	Product Description	Brand	Inner Pack	Outer Pack
1125	6ml Pre-filled Syringe	OptiLube™	Pack of 25	Case of 250
1126	11ml Pre-filled Syringe	OptiLube™	Pack of 25	Case of 250

Syringes are designed specifically for precise application. Each pre-filled syringe is sterile and individually wrapped ensuring minimal handling, maintaining an aseptic and 'gold-standard' procedure. OptiLube™ sterile lubricating jelly is available in 6ml and 11ml pre-filled syringes.

- Measuring facility on side of syringe enables precise application of a controlled volume of jelly
- Latex free bung and plunger
- Easy to prime syringe
- Manufacturing, expiry and LOT numbers printed on each individual syringe's tear off strip. These are easily retained for accurate recording of patients' records



Tubes 5g e 42g e 82g e 113g e

OMS Code	Product Description	Brand	Inner Pack	Outer Pack
1124	5g Tube	OptiLube™	Pack of 48	Case of 480
1121	42g Tube	OptiLube™	Pack of 12	Case of 72
1122	82g Tube	OptiLube™	Pack of 12	Case of 72
1127	113g Tube	OptiLube™	Pack of 12	Case of 72

The classic choice for clinicians. The tubes have a special design flip-cap. This feature makes the tubes easy to open and avoids spillage. Available in 5g mini tubes, 42g, 82g and a large 113g tube.

- Pliable plastic coated tubes avoids wastage as clinician can extract as much gel as possible
- Large cap enables the tubes to be stored 'standing up' on the cap reducing storage space



"Sheffield Teaching Hospitals NHS Foundation Trust changed to Optilube sterile 5g sachet back in November 2010. The Clinicians like the small individual sachets. The gel has the right consistency and suits all our clinical needs. We are very happy with the product."

Michaela Fairest, Clinical Specialist,
Sheffield Teaching Hospitals NHS Trust

"We switched our whole group of hospitals to OptiLube."

Angela Timms, Procurement Team,
United Lincolnshire Hospitals NHS Trust

"Choosing the OptiLube range has been a great decision for my practice and patients. The gel is high quality and consistently performs extremely well on all of the many procedures it is required for."

Dr Donna Kershner, Doctor, York PCT, UK

ENVIRONMENT POLICY



As a manufacturer of medical products we recognise that we have a responsibility to reduce our impact on the environment and to conserve natural resource. We endeavour to use the smallest viable packaging size to minimise material usage and therefore landfill footprint.



When transporting products we make sure our transport operations are as efficient and environmentally friendly as possible.

All printed marketing materials are 50% recycled and certified under FSC's Chain-of-Custody programme.



Office: 47 Park Square East, Leeds, LS1 2NL, UK
Head Office: Malmarc House, 116 Dewsbury Road, Leeds, LS11 6XD
Tel.: +44 (0) 8456 435 479, Fax: +44 (0) 1142 383 826
www.optimummedical.co.uk