

# SafeGuide® Over the Guidewire Esophageal Dilators



## Flexibility and control for exacting procedures.

At times, there is need for greater accuracy. That's why there's SafeGuide®, a remarkably flexible over the guidewire solution that ensures exact placement of the dilator. Position the guidewire with the assistance of an endoscope and direct the SafeGuide® dilator where you need it.

### Benefits

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Latex-free
- Dilators are reusable

### Features

- 16 sizes: 15 FR to 60 FR (in 3 FR increments)
- Choose a full set, mini set or individual pieces
- Easy to clean with optional brushes and adapters
- With proper care, reusable SafeGuide® dilators provide years of service



Store and transport SafeGuide® dilators in a convenient storage case. A labeled slot holds each dilator secure. The case holds a full set of 16 dilators.

### SafeGuide® Dilatation Systems

Set	Part No.	Contents
16-piece full set	1214-00	16 dilators, 15-60 FR, 1 storage case, 1 marked spring tip guidewire, 1 reusable cleaning brush, 25 disposable cleaning adapters
7-piece mini set	1214-10	7 dilators (15, 21, 27, 33, 39, 42, 45 FR), 1 storage case, 1 marked spring tip guidewire, 1 reusable cleaning brush, 25 disposable cleaning adapters
Custom mini set	Contact your sales representative	7 dilators of your choice, 1 storage case, 1 marked spring tip guidewire, 1 reusable cleaning brush, 25 disposable cleaning adapters



The single-piece construction of the SafeGuide® Spring Tip Guidewire resists breakage and flexes easily.



### SafeGuide® Over the Guidewire Esophageal Dilators

FR (French) Size	(mm)	Part No.
15	5	1214-15
18	6	1214-18
21	7	1214-21
24	8	1214-24
27	9	1214-27
30	10	1214-30
33	11	1214-33
36	12	1214-36
39	13	1214-39
42	14	1214-42
45	15	1214-45
48	16	1214-48
51	17	1214-51
54	18	1214-54
57	19	1214-57
60	20	1214-60


### Accessories



Storage case indicates size.

### SafeGuide® Accessories

Set	Part No.	Description	Qty/Case
SafeGuide® Storage Case	1214-01	Holds up to 16 SafeGuide® Dilators. 36"L x 10"W x 4"H	1
Marked Spring Tip Guidewires	1214-02	Guides the dilator. 210cm	2
Disposable Cleaning Brushes	1214-03D	Cleans the inner lumen of the dilator	3
Reusable Cleaning Brushes	1214-03R	Cleans the inner lumen of the dilator	3
Disposable Cleaning Adapters	1214-04	Used to flush the inner lumen of the dilator	25



Clearly marked reusable and single-use cleaning brushes.



Insert cleaning brush into the distal end and push through to the other side for complete inner lumen cleaning.



Disposable cleaning adapter.

