

HS TRAPSYSTEM® SET

BONE MARROW BIOPSY SET WITH TRAPSYSTEM DEVICE

The **TRAPSYSTEM®** grants a quality marrow sampling by trapping the specimen into the needle cannula without damaging it and without performing the painful needle deflection technique.

- Greater patient comfort
- Safer marrow sampling
- Simpler manoeuvrability
- Faster procedure

The needle, with its concave cutting edge point, allows to easily penetrate into the bone marrow cavity. Also the marked stylet, with its pyramidal point, allows to easily penetrate into the bone marrow cavity. The ergonomically designed handle enables a safer manoeuvrability, since the strength needed to penetrate the bone marrow cavity is homogeneously distributed on the whole palm of your hand and not locally.

CODE	Diameter	Length	Box
TRAPJ 0810	8G	100mm	10
TRAPJ 0815	8G	150mm	10
TRAPJ 1110	11G	100mm	10
TRAPJ 1115	11G	150mm	10
TRAPJ 1306	13G	65mm	10
TRAPJ 1310	13G	100mm	10



H.S. Hospital Service S.p.A. is certified ISO 9001 and EN 46001

BONE MARROW BIOPSY



H.S. HOSPITAL SERVICE S.p.A.

Via delle Valli snc - 04011 Aprilia (LT)
Tel. 069201961 - Fax 0692727871 - Orders fax 069275519
www.hshospitalservice.com e-mail: hs@hshospitalservice.com