


# Bougie Dilator Sets



CE 0123

 常美医疗  
Changmei Medtech

Concentration Makes Profession

# Bougie Dilator Sets

## Usage

- It is used for the expansion and treatment of esophagus and cardia.

## Characteristics

- Made by medical polymer materials with proper softness.
- Anti-shedding design, more safety and reliable.
- The front cone angle is manufactured by mechanized processing. With smooth transition, it can be inserted into the narrow parts successfully.
- Two radiopaque marks on both ends make it easy to observe under X-ray.
- Print scales on both sides, allowing the physicians to identify the insertion length.

## Specifications(Unit: mm)

Model	Dilator Diameter	Working Length	Guide Wire
BDS-05	5	850	0.035"
BDS-06	6	850	0.035"
BDS-07	7	850	0.035"
BDS-08	8	850	0.035"
BDS-09	9	850	0.035"
BDS-10	10	850	0.035"
BDS-11	11	850	0.035"
BDS-12	12	850	0.035"
BDS-13	13	850	0.035"
BDS-14	14	850	0.035"
BDS-15	15	850	0.035"
BDS-16	16	850	0.035"
BDS-17	17	850	0.035"
BDS-18	18	850	0.035"
BDS-19	19	850	0.035"
BDS-20	20	850	0.035"



Tel: 0086-519-88520118

Fax: 0086-519-88551179-8888

Web: [www.czmed.com](http://www.czmed.com)

E-mail: [sales@czmed.com](mailto:sales@czmed.com)

Add: No. 27 Xinke West Road, Luoyang Town, Wujin District, Changzhou, Jiangsu, China

