



Interventional Radiology

www.allwinmedical.com

INDEX

INFORMATION	1
DRAINAGES	2
GUIDEWIRES	3
DILATORS	4
NEEDLES	5
BIOPSY NEEDLES	6
SPECIAL RADIOLOGY	7

PRODUCT WARRANTY

1

Allwin warrants that at the time of manufacture, the products are prepared and tested in accordance with highest standards. Due to biological differences in individuals, no product is 100% effective under all circumstances for all patients. Further, because we have no control of the condition under which the product is used, diagnosis of patient the method of use or administration, and handling and storage of the product after it leaves our possession, we do not warrant any consequence which may arise pertaining to the final use of the device. The foregoing warranty is exclusive and in lieu of all other warranties written, oral or implied. No representative of Allwin may change any of the foregoing and the buyer accepts the product subject to all terms hereof.

DESIGN CHANGES

Continuous improvements is the chosen path at Allwin and we reserve the right to change or discontinue any product without notice.

RETURN GOODS POLICY

Allwin will accept returned goods only after issuing a Return Goods Authorisation (RGA) to the requester within 30 days of purchase. The request for return must be accompanied by reason for return invoice number, order number. Credit will be issued as per company rules effective at the time of return and mentioned in RGA. Special, opened & damaged products are not returnable.

SPECIAL ORDERS

Allwin welcomes new ideas and will manufacture new or special items, if within the capability of the company. Modification to the original design of any product can be made with only a slight delay in delivery schedule


STORAGE OF PRODUCTS

Allwin recommends that the product be stored in cool and dry place. Plastic material used may degrade following extended periods of exposure to light and temperature.



Ultima Drainage Sets

Used for percutaneous placement of a pigtail catheter in the pelvis for nephrostomy drainage.



2.1

- A convenient complete set for percutaneous placement of drainage catheter - ultima pigtail design.
- Material provides superior biocompatibility and drainage.
- Highly radiopaque for precise placement.
- Locking loop design available for better anchoring.

Set Includes :

- Initial Puncture Needles - 3 Part. • J Tip, Fixed Core, Guidewire • Fascial Dilators • Ultima Pigtail Drainage Catheter • Urine Bag Connector • Scalpel • Catheter fixing disc • Stop cock

QTY/Box: 1

Order No. 30 cm	Order No. 20 cm	With Fascial Dilators	Fr
U-PDS-P8/30	U-PDS-P8/20	7,8,9	8
U-PDS-P9/30	U-PDS-P9/20	7,8,9,10	9
U-PDS-P10/30	U-PDS-P10/20	6,8,10,12	10
U-PDS-P12/30	U-PDS-P12/20	6,8,10,12,14	12
U-PDS-P14/30	U-PDS-P14/20	6,8,10,12,14,16	14

Combinations:

With Locking Loop System: Add Suffix - LS

With Lunderquist guidewire instead of PTFE guidewire: Add Suffix - L

With Hydrophilic coating: Add Suffix - H


e.g.

U-PDS-P7-LS

U-PDS-P7-L

Ultima Drainage Catheter

2.2



- A convenient complete set for percutaneous placement of drainage catheter - ultima pigtail design.
- Material provides superior biocompatibility and drainage.
- Highly radiopaque for precise placement.
- Locking loop design available for better anchoring.

QTY/Box: 1

Order No. 30 cm with UBC	Order No. 20 cm with UBC	Order No. 30 cm	Order No. 20 cm	Fr
U-PCN-P8/30-U	U-PCN-P8/20-U	U-PCN-P8/30	U-PCN-P8/20	8
U-PCN-P9/30-U	U-PCN-P9/20-U	U-PCN-P9/30	U-PCN-P9/20	9
U-PCN-P10/30-U	U-PCN-P10/20-U	U-PCN-P10/30	U-PCN-P10/20	10
U-PCN-P12/30-U	U-PCN-P12/20-U	U-PCN-P12/30	U-PCN-P12/20	12
U-PCN-P14/30-U	U-PCN-P14/20-U	U-PCN-P14/30	U-PCN-P14/20	14

Combinations:

With Locking Loop System: Add Suffix - LS

With Hydrophilic coating: Add Suffix - H

Ultima Direct Puncture Drainage Set

Used for percutaneous nephrostomy following the "direct-puncture" technique particularly when there is a dilated urine collection system, a hydronephrosis.



2.3

- The catheter is back-loaded on a two-part Trocar Needle and straightened and then introduced through a skin incision directly.
- The procedure does away with the use of guidewire and fascial dilators as there is no dilatation performed.
- The set may also be used for cyst and abscess drainage.
- With Stop cock.

QTY/Box: 1

Order No. 30 cm	Order No. 20 cm	Fr
U-DPD-8/30	U-DPD-8/20	8
U-DPD-9/30	U-DPD-9/20	9
U-DPD-10/30	U-DPD-10/20	10
U-DPD-12/30	U-DPD-12/20	12
U-DPD-14/30	U-DPD-14/20	14

Combinations:

With Locking Loop System: Add Suffix - LS

With Hydrophilic coating: Add Suffix - H

DRAINAGES

Percutaneous Nephrostomy Drainage Set - Pigtail

Used for percutaneous placement of a pigtail catheter in the pelvis for nephrostomy drainage.

2.4



- A convenient complete set for percutaneous placement of drainage catheter - pigtail design.
- Material provides superior biocompatibility and drainage.
- Highly radiopaque for precise placement.
- Locking loop design available for better anchoring.
- Softer polyurethane

Set Includes :

- Initial Puncture Needles - 3 Part . • J Tip, Fixed Core, PTFE Coated 0.038", Guidewire • Fascial Dilators
- Pigtail Drainage Catheter 30cm • Urine Bag Connector • Scalpel • Catheter fixing disc • Stop cock

QTY/Box: 1

Order No.	With Fascial Dilators	Fr
PDS-P7	6, 7, 8	7
PDS-P8	7, 8, 9	8
PDS-P9	7, 8, 9, 10	9
PDS-P10	6, 8, 10, 12	10

Combinations:

With Locking Loop System: Add Suffix - LS

With Lunderquist guidewire instead of PTFE guidewire: Add Suffix - L

With Hydrophilic coating: Add Suffix - H


e.g.

PDS-P7-LS

PDS-P7-L

Percutaneous Nephrostomy Drainage Set - Malecot

Used for percutaneous placement of a malecot catheter in renal pelvis for nephrostomy drainage.



2.5

- A convenient complete set for percutaneous placement of drainage catheter - malecot design. -ideally suitable for the drainage of thick, viscous fluids.
- Material provides superior biocompatibility and drainage.
- Highly radiopaque for precise placement.
- 30 cm.

Set includes :

- Initial Puncture Needles - 3 Part • J Tip, Fixed Core, PTFE Coated 0.038", 150 cms. Guidewire • Fascial Dilators • Malecot Drainage Catheter • Urine Bag Connector • Scalpel • Catheter fixing disc • Stop cock

QTY/Box: 1

Order No.	With Fiscial Dilators	Fr
PDS-M10	6, 8, 10, 12	10
PDS-M12	6, 8, 10, 12, 14	12
PDS-M14	6, 8, 10, 12, 14, 16	14
PDS-M16	6, 8, 10, 12, 14, 16, 18	16

Combinations:

With Lunderquist guidewire instead of PTFE guidewire: Add Suffix - L


With Hydrophilic coating: Add Suffix - H

e. g.

PDS-M10-L

Percutaneous Nephrostomy Catheter - Pigtail

2.6



- Pigtail design.
- Biocompatible, radiopaque catheter.
- With Stop cock.
- Contact for other sizes.

QTY/Box: 10

Order No. with Urine Bag Connector	Order No.	Length (cm)	Fr
PCN-P7-U	PCN-P7	30	7
PCN-P8-U	PCN-P8		8
PCN-P9-U	PCN-P9		9
PCN-P10-U	PCN-P10		10


Combinations:

With Locking Loop System: Add Suffix - LS

e.g.

PCN-P7-LS

Percutaneous Nephrostomy Catheter - Malecot




- Malecot design ideally suited for thick, viscoous fluids.
- Biocompatible, radiopaque catheter.
- With flexible stylet.
- With Stop cock.

QTY/Box: 10

Order No. with Urine Bag Connector	Order No.	Length (cm)	Fr
PCN-M10-U	PCN-M10	30	10
PCN-M12-U	PCN-M12		12
PCN-M14-U	PCN-M14		14
PCN-M16-U	PCN-M16		16

SupraSet DPD

Used for percutaneous bladder drainage.



• 30 cm.

2.7

Set includes :

• Drainage Catheter • Needle • Urine Bag Connector • Stop Cock


QTY/Box : 5

Order No.	Fr
SS-DPD 12	12
SS-DPD 14	14
SS-DPD 16	16

Perc Drain Catheter Set

Used for nephrostomy drainage through an existing track

2.8




Set includes :

- Drainage Catheter
- Stylet
- Urine Bag Connector
- Catheter Fixing disc
- Stop cock

QTY/Box :5

Order No.	Length	Fr
PDCS 18	30	18
PDCS 22		22
PDCS 24		24

LADS - Liver Abscess Drainage Set



- A convenient complete set for percutaneous placement of drainage catheter - malecot design ideally suitable for the drainage of thick viscous fluids.
- Material provides superior biocompatibility and drainage.
- Highly radiopaque for precise placement.

2.9

QTY/Box: 1

Order No.	Fr	With Fascial Dilators
LADS – M10	10	6, 8, 10, 12
LADS – M12	12	6, 8, 10, 12, 14
LADS – M14	14	6, 8, 10, 12, 14, 16
LADS – M16	16	6, 8, 10, 12, 14, 16, 18

Set Includes: • Initial Puncture Needle - 3 Part • J Tip Fixed Core, PTFE Coated 0.038" 150 cm Guidewire • Fascial Dilators • Malecot Drainage Catheter with Stylet • Urine Bag Connector • Scalpel


Combinations:

With Lunderquist guidewire instead of PTFE guidewire: Add Suffix - L

Re-Entry Malecot Catheter

Used for internal as well as external drainage with combination of nephrostomy catheter with stent.

2.10



- Material provides superior biocompatibility and drainage.
- Highly radiopacity.
- Tip design allows guidewire to be passed in ureter.


QTY/Box: 1

Order No with Urine Bag Connector	Order No.	GW	Length (cm)	Fr
RM 10-U	RM 10	0.038"	30	10
RM 12-U	RM 12			12
RM 14-U	RM 14			14
RM 16-U	RM 16			16
RM 18-U	RM 18			18
RM 20-U	RM 20			20

J Drain Percutaneous Drainage Set

Used for abscess and other multipurpose drainage using seldinger technique.

2.11



- A convenient complete set for percutaneous placement of drainage catheter - J drainage design.
- Material provides superior biocompatibility and drainage.
- Highly radiopaque for precise placement.
- Locking loop design available for better anchoring.

Set Includes :

- Initial Puncture Needles - 3 Part. • J Tip, Fixed Core, Guidewire • Fascial Dilators • Ultima Pigtail Drainage Catheter • Urine Bag Connector • Scalpel • Catheter fixing disc • Stop cock

QTY/Box : 1

Order No. 30 cm	Order No. 20 cm	With Fascial Dilators	Fr
J-PDS-8/30	J-PDS-8/20	7,8,9	8
J-PDS-9/30	J-PDS-9/20	7,8,9,10	9
J-PDS-10/30	J-PDS-10/20	6,8,10,12	10
J-PDS-12/30	J-PDS-12/20	6,8,10,12,14	12
J-PDS-14/30	J-PDS-14/20	6,8,10,12,14,16	14

Combinations:

With Locking Loop System: Add Suffix - LS

With Lunderquist guidewire instead of PTFE guidewire: Add Suffix - L

With Hydrophilic coating: Add Suffix - H

e.g.

U-PDS-P7-LS

U-PDS-P7-L

J Drain Direct Puncture Drainage Set

Used for abscess and other multipurpose drainage using seldinger technique.

2.12



- The procedure does away with the use of guidewire and fascial dilators as there is no dilatation performed.
- The set may also be used for cyst and abscess drainage.
- With Stop cock.
- 30 cm.

QTY/Box : 1

Order No.	Fr
J-DPD-8	8
J-DPD-9	9
J-DPD-10	10
J-DPD-12	12
J-DPD-14	14


Combinations:

With Lunderquist guidewire instead of PTFE guidewire: Add Suffix - L

With Hydrophilic coating: Add Suffix - H

PTCD - Percutaneous Transhepatic Cholangiography Drainage

Used to perform cholangiography and for internal / external bile duct drainage.



- Complete set for percutaneous drainage.
- Biocompatible, radiopaque catheter.
- Available with Straight and Special pigtail.

2.13

QTY/Box: 5

Order No.	Size (Fr)
PTCD – S8	8
PTCD – S9	9
PTCD – S-K8	8
PTCD – S-K9	9

QTY/Box: 5

Order No.	Size (Fr)
PTCD – SP8	8
PTCD – SP9	9
PTCD – SP-K8	8
PTCD – SP-K9	9


Set includes :

- Straight / Special Pigtail Catheter
 - Initial Puncture Needle - 3 Part
 - Fascial Dilators
 - Drainage Connector
 - Scalpel
 - Catheter Fixing Disc
 - Guidewire - PTFE
- * With 8Fr Kit, 7Fr, 8Fr, 9Fr, Fascial Dilators available
 * With 9Fr Kit, 7Fr, 8Fr, 9Fr, 10Fr Fascial Dilators available
 * With Lunderquist guidewire instead of PTFE guidewire: Add Suffix - L

Centesis Catheter Needle

Used for percutaneous fluid aspirations.

2.14




- With 4 side holes
- Radiopaque catheter

QTY/Box: 5

Order No.	Fr	Length (cm)	Needle Gauge
CCN 4/7	4	7	20
CCN 4/10	4	10	20
CCN 5/7	5	7	19
CCN 5/10	5	10	19
CCN 5/15	5	15	19

HyGuide Hydrophilic Guidewires

For percutaneous entry into vessel using the seldinge technique or for endoscopic access.



- Exacting tolerances of each wire facilitates easy reliable instrument delivery.
- Stainless steel shaft with hydrophilic coating.


3.1

QTY/Box: 5

Order No.	Description
HGSF 150-35	Straight fixed core, Hydrophilic coated 150cm-0.035"

J - Tip Fixed Core Teflon Guidewire


3.2



QTY/Box: 5

Order No.	Description
JFT 150 - 32 - 3	J tip, 3mm, Fixed Core, Teflon Coated, 150cm 0.032"
JFT 150 - 35 - 3	J tip, 3mm, Fixed Core, Teflon Coated, 150cm 0.035"
JFT 150 - 38 - 3	J tip, 3mm, Fixed Core, Teflon Coated, 150cm 0.038"

Straight Tip Fixed Core Teflon Guidewire



QTY/Box: 5

Order No.	Description
SFT 150-25	Straight, Fixed Core, Teflon Coated, 150cm 0.025"
SFT 150-28	Straight, Fixed Core, Teflon Coated, 150cm 0.028"
SFT 150-32	Straight, Fixed Core, Teflon Coated, 150cm 0.032"
SFT 150-35	Straight, Fixed Core, Teflon Coated, 150cm 0.035"
SFT 150-38	Straight, Fixed Core, Teflon Coated, 150cm 0.038"
SFT 400-35	Straight, Fixed Core, Teflon Coated, 400cm 0.035"
SFT 500-35	Straight, Fixed Core, Teflon Coated, 500cm 0.035"

Schuller Guidewires

For percutaneous entry into vessel using the seldinger technique or for endoscopic access.



- For easy catheter manipulation and introduction of instruments.
- Catheters and other instruments can easily be introduced over the floppy part connected to the rigid shaft.

3.3

QTY/Box: 1

Order No.	Shaft	Length (cm)	Rigid Shaft	Floppy Part	Tip
SHJ 100-30-60	.035"	100	30cm	60 cm	J-curved, 3.0mm

Lunderquist Guidewires

For percutaneous entry into vessel using the seldinger technique or for endoscopic access.



- Firm and easy to introduce
- Shapeable tip facilitates accurate placement even in difficult anatomy.
- Available in 100 cm length
- 3 mm J-Tip

QTY/Box: 1


Order No.	Length (cm)	Diameter tip	Tip
LQJ-100-35-3	100	0.035"	3mmJ

DILATORS

Amplatz Renal Dilator Set

Amplatz Renal Dilator is used for progressive dilatation for the nephrostomy tract prior to percutaneous kidney stone removal.

4.1



- Radiopaque dilators & sheaths.
- Short tapered tips.
- Smooth surface reduces patient's trauma.
- Introducer Catheter acts as a guide for dilators 14F - 30F.
- Introducer Catheter Sheath allows safety wire placement.
- 30 cm.

Amplatz Full Set

QTY/Box: 1

Order No.	Description
ARDS 8/30 – FS	<p>Amplatz renal Dilatation Full Set Includes:</p> <ol style="list-style-type: none"> 1) 8, 10, 12, 14, 16, 18 Fr – PE Facial Dilators 2) 20,24, 26, 28, 30 Fr – PE Renal Dilators 3) 24, 26, 28, 30 Fr Amplatz Sheath – PTFE 4) Introducer Catheter for Dilators 8 Fr, 70cm 5) Introducer Catheter Sheath PTFE

Amplatz Set

QTY/Box: 1

Order No.	Description
ARDS 8/30	<p>Amplatz renal Dilatation Set Includes:</p> <ol style="list-style-type: none"> 1) 8, 10, 12, 14, 16, 18 Fr – PE Facial Dilators 2) 20,24, 26, 28, 30 Fr – PE Renal Dilators 3) 28, 30 Fr Amplatz Sheath – PTFE 4) Introducer Catheter for Dilators 8 Fr, 70cm

DILATORS

Amplatz Renal Dilator




QTY/Box: 5

Order No.	Size (Fr)	Length (cm)	Description
RD 20	20	30	Renal Dilator 20 F
RD 22	22		Renal Dilator 22 F
RD 24	24		Renal Dilator 24 F
RD 26	26		Renal Dilator 26 F
RD 28	28		Renal Dilator 28 F
RD 30	30		Renal Dilator 30 F

4.2

Renal Dilator with Sheath



QTY/Box: 5

Order No.	Size (Fr)	Length (cm)	Description
RD-AM 20/30	20	30	Renal Dilator with Sheath 20Fr, 30cm
RD-AM 24/30	24		Renal Dilator with Sheath 24Fr, 30cm
RD-AM 26/30	26		Renal Dilator with Sheath 26Fr, 30cm
RD-AM 28/30	28		Renal Dilator with Sheath 28Fr, 30cm
RD-AM 30/30	30		Renal Dilator with Sheath 30Fr, 30cm

DILATORS

Amplatz Introducer Sheath




4.3

QTY/Box: 5

Order No.	Size (Fr)	Length (cm)	Description
ICS 10/30	10	30	Introducer Catheter Sheath PTFE

Introducer Catheter




QTY/Box: 5

Order No.	Size (Fr)	Length (cm)	Description
IC 8/70	8	70	Introducer Catheter

DILATORS

Fascial Dilators

Fascial dilators are used for dilatation of percutaneous track over the guidewire.



- Fully radiopaque polyethylene.
- Female luer lock allows heparin, contrast media and other solutions to be injected using standard luer lock syringe.
- Accepts 0.038" guidewires.


4.4

QTY/Box: 5

Order No.	Length (cm)	Description
FD 6	20	Fascial Dilator 6F
FD 7		Fascial Dilator 7F
FD 8		Fascial Dilator 8F
FD 9		Fascial Dilator 9F
FD 10		Fascial Dilator 10F
FD 12		Fascial Dilator 12F
FD 14		Fascial Dilator 14F
FD 16		Fascial Dilator 16F

Safety Wire Introducer

Used for placement of safety wire into the pelvis.



QTY/Box: 10


Order No.	Length (cm)	Fr
SWI	20	10/6

DILATORS

Biliary Dilatation Catheter

Used to dilate the biliary strictures, pancreatic strictures or papillary orifice.

4.5



• Length 200 cm.

BDC [PEBAX]

QTY/Box: 5

Order No.	Transition Length	TAPER	
		From French Size	To French Size
BDC 7/5	2.5 cm	5F	7F
BDC 10/6	2.5 cm	6F	10F

BDC-T [PTFE]

QTY/Box: 5

Order No.	Transition Length	TAPER	
		From French Size	To French Size
BDC 7/5 -T	2.5 cm	5F	7F

NEEDLES

Chiba Needle -2- Part Bevel Tip

Used for aspiration and cytological biopsy.



5.1

QTY/Box: 5

Order No.	G	Length (cm)	Tip
CH-22-15-2	22	15	2
CH-22-20-2	22	20	2
CH-22-28-2	22	28	2

NEEDLES

PTC Needle



- With sheath
- Trocar point

5.2

QTY/Box: 5

Order No.	Fr	Length (cm)
PTCN 5/15	5	15
PTCN 5/25	5	25

Initial Puncture Needle -3- Part Bevel Tip

Used for percutaneous nephrostomy and all other percutaneous application where a 0.035 " guidewire is to be used.



QTY/Box: 5

5.3

Order No.	Description
IPN-18-20-3	Initial Puncture Needle, 18G, 20cm, 3 part Bevel

Initial Puncture Needle -2- Part Trocar Tip

Used for percutaneous nephrostomy and all other percutaneous application where a 0.035" guidewire is to be used.



QTY/Box: 5

Order No.	G	Length (cm)	Tip
IPN-18-20-2T	18	20	2T
IPN-18-30-2T	18	30	2T

NEEDLES

TLA Needle

Used for placement of 0.038" guidewire.




- Radiopaque sheath for better visualization.
- Trocar point for easy introduction
- 5 Fr sheath

5.4

QTY/Box: 5

Order No.	G	Length (cm)
TLA-19	19	20

Introducer Needle




QTY/Box: 5

Order No.	G	Length (cm)
IN-18/7	18	7

Breast Localization Needle

Used for pre-surgery localization of breast non-palpable lesions.



- Easy thread wire entry.
- Accurate depth placement with centimeter marks.
- Wire is marked in two places: the mark at the leading edge indicates the wire is at the distal tip of needle and the mark at the end indicates the wire deployment.
- Needle hub design allows dye injection or aspiration of fluid.


5.5

QTY/Box: 5

Order No.	G	Length (cm)
BLN 20/4	20	4
BLN 20/7	20	7
BLN 20/9	20	9

Speedo Biopsy Needle

Used for soft tissue biopsy.



- User friendly, lightweight, compact and provided with a strong mechanism: it is therefore more reliable in cutting the most compact tissue.
- Provided with 2 triggers; one in the front and one in the back to facilitates operations in every position.
- Safety: the safety switch is automatically activated when loading the gun.
- Centimeter marking facilitate precise depth placement.
- The tip is treated to make it echogenic and facilitate its positioning under ultrasound guidance.
- Available at select places only

6.1

QTY/Box: 5

Order No.	Length (cm)	Gauge	Description
SP 18/16	16	18	Speedo, 18G, 16cm.
SP 18/20	20	18	Speedo, 18G, 20cm.
SP 18/25	25	18	Speedo, 18G, 25cm.

BIOPSY NEEDLES

Ultracore Biopsy Needle - Manual

Ultracore is used for performing core biopsies of soft tissues.



- Ultrasonically visible tip facilitates precise placement.
- CM marking for proper positioning.
- Wide guage and length range.

QTY/Box: 5

6.2

Order No.	G	Length (cm)
UCB 14/3 UCB 14/4.5 UCB 14/6	14	3" 4.5" 6"
UCB 16/3 UCB 16/4.5 UCB 16/6	16	3" 4.5" 6"
UCB 18/3 UCB 18/4.5 UCB 18/6	18	3" 4.5" 6"

Bone Marrow Biopsy Needle



- Stylet with trocar edge fast and easy penetration.
- Ergonomically designed handle enables a comfortable procedures.
- The adjustable and removable needle stop, allows needle length variation with safety controlling the penetration depth.


6.3

QTY/Box: 5

Order No.	Length (cm)	Gauge (G)
BMN 11/10	10	11
BMN 13/7	7	13

Hystero SonoSalpinography Catheter

HSG Catheter is used for injection of contrast media into the uterine cavity to evaluate patency of the fallopian tubes.



- Strong balloon with 2 Stop cock & 1 syringe
- Closed end

QTY/Box: 1

Order No.	Description
HSG	Hystero SonoSalpingography Catheter


7.1

Caution:

This product contains natural rubber latex which may cause allergic reaction in some individuals.

Fallopian Tube Recanalisation Set - FTR-M

Used for fallopian tube recanalisation under radiological fluoroscopic control.



- 5 Fr Pebax Catheter - preshaped
- Malleable stylet
- Guidewire

QTY/Box : 1

Order No.	Description
FTR - M	Fallopian Tube Recanalisation Set

7.2

Fallopian Tube Recanalisation Set - FTR

Used for fallopian tube recanalisation under radiological fluoroscopic control.




- 5 Fr selective injection catheter
- 3 Fr inner catheter enters tubal OS
- 0.018" guidewire

QTY/Box : 1

Order No.	Description
FTR	Fallopian Tube Recanalisation Set

Small Bowel Enema Catheter Open / Closed Tip

Used for double contrast display of small intestine.



- Round tip for atraumatic introduction.
- Soft radiopaque, flexible catheter for quick introduction.
- Heavy duty torque wire with manipulation knob.

Small Bowel Enema Catheter Closed Tip

QTY/Box: 1

Order No.	Description
SBEC-C 12/6-G	SBEC, Closed Tip, 12F, 6 sides holes with Torque wire, 150 cm
SBEC-C 14/6-G	SBEC, Closed Tip, 14F, 6 sides holes with Torque wire, 150 cm
SBEC-C 12/6	SBEC, Closed Tip, 12F, 6 sides holes 150 cm
SBEC-C 14/6	SBEC, Closed Tip, 14F, 6 sides holes 150 cm

7.3

Small Bowel Enema Catheter Open Tip

QTY/Box: 1

Order No.	Description
SBEC-O 12/0-G	SBEC, Open Tip, 12F, 0 sides holes with Torque wire, 150 cm
SBEC-O 14/0-G	SBEC, Open Tip, 14F, 0 sides holes with Torque wire, 150 cm
SBEC-O 12/0	SBEC, Open Tip, 12F, 0 sides holes 150 cm
SBEC-O 14/0	SBEC, Open Tip, 14F, 0 sides holes 150 cm

Inflation Device Advanced



- Ergonomic design facilitates excellent eye to hand co-ordination.
- Efficient Intra cam locking mechanism allows to maintain pressure and guarantee fast and automatic release.
- Superior design with grained surface and slip free lever ensures maximum grip even with wet hands.
- Pressure gauge ensures fatigue free readability.
- Up to 30 / 26 ATM Analogue pressure range.
- Registers immediately without pressure decay or drift.

QTY/Box: 1

Order No.	Description
IDA 20	Inflation Device Advanced 20 cc

7.4

ACE Inflation Device



- Ergonomic design facilitates excellent eye to hand co-ordination.
- Efficient Intra cam locking mechanism allows to maintain pressure and guarantee fast and automatic release.
- Superior design with grained surface and slip free lever ensures maximum grip even with wet hands.
- Pressure gauge ensures fatigue free readability.
- Up to 30 / 26 ATM Analogue pressure range.
- Registers immediately without pressure decay or drift.

QTY/Box: 1

Order No.	Description
AID 20	Ace Inflation Device 20cc

Drainage Bag Connector



QTY/Box: 25

Order No.	Description
DBC	Drainage Bag Connector
DBC-SC	Drainage Bag Connector with stop cock
DBC-S	Drainage Bag Connector short

7.5

One Way Stopcock



QTY/Box: 50

Order No.	Description
OWS	One Way Stopcock

Catheter Fixing Disk




- Used to stabilize indwelling catheter.
- Secured loading with pull-tie.
- Two Sizes.

QTY/Box: 25

Order No.	Description
CFD 12	Upto 12F catheter
CFD 20	Upto 20F catheter

7.6

Scalpel



QTY/Box: 50

Order No.	Description
SP	Scalpel

