



The j carina stent is used to bridge an anastomosis leakage after pneumonectomy. Because of the matter of urgency we offer you different sizes of the j carina stents.

The j carina stent is preloaded on the reload delivery system. For correction of the stent position the stent can be reloaded into the delivery system during the procedure. The excellent tantalum marker makes the stents clearly visible in the fluoroscopy.

- Stent with complete silicone covering
- Atraumatic, rounded ends
- Memory effect of nitinol material – high radial force and flexibility
- Radiopaque tantalum marker
- Reload delivery system

Length 600 mm, OTW  
 Ø 8 mm / 24 French  
 Accepts guide wire 0.035 inch

aerstent · j carina · complete silicone covered

REF	Trachea Ø	Trachea Length	Bronchus Ø	Bronchus Length
533-16-030	16	30	12	30
533-18-040	18	40	14	30
533-20-050	20	50	14	30