


E.R.C.P. Catheters

Single Use

E.R.C.P. Catheters


- "Pre-Curved" Catheter with Memory Effect
Easier Positioning to the Papilla and Possible Manual Tip Manipulation for Individual Anatomy
- Guide Wire Compatibility
Contrast Injection while Guide Wire is in Place
- Metal Tips / Platinum Rings
Maximum X-Ray Visualization
- Color Graduation
Macroscopically Controlled Cannulation Depth



E.R.C.P. Catheters

Single Use

E.R.C.P. Catheters



Compatible to .035" Guide Wire

E.R.C.P. Catheter – “Pre-Curved” with Color Graduation and Metal Tips / Platinum Rings, Removable Nitinol Stylet, sterile packed

Catalog No.	Tip Shape	Catheter Outer Diameter
GER-01-18-035	Standard Tip	1.8 mm
GER-02-18-035	5 mm Tapered Tip	1.8 mm
GER-04-18-035	Metal Cone Tip	1.8 mm
GER-05-18-035	Round Metal Tip	1.8 mm
GER-06-18-035	Bullet Type Metal Tip	1.8 mm
GER-07-18-035	Metal Ball Tip	1.8 mm

Compatible to .021" Guide Wire

E.R.C.P. Catheter – “Pre-Curved” with Color Graduation and Platinum Ring, Removable Nitinol Stylet, sterile packed

Catalog No.	Tip Shape	Catheter Outer Diameter
GER-03-18-021	10 mm Tapered Tip	1.8 mm

Distributed by: