


Nasal Biliary Drainages


Single Use

Individual Components


Nasal Tube

Catalog No.	French Size	Length	Product
GND-41-05-050	5.0	50 cm	
GND-41-06-050	6.0	50 cm	
GND-41-07-050	7.0	50 cm	
GND-41-10-050	10.0	50 cm	

Multi-Medfix Injection Connector

Catalog No.	French Size	Product
GND-43-05-000	5.0	
GND-43-06-000	6.0	
GND-43-07-000	7.0	
GND-43-10-050	10.0	

Connecting Tube

Catalog No.	Length	Product
GND-42-03-135	135 cm	

Distributed by:

Cobra Medical BV
 Stavangerweg 21-14, 9723 JC Groningen
 Postbus 5153, 9700 GD Groningen
 Nederland
 T +31 (0)50 52 80 498
 F +31 (0)50 52 80 495



www.cobramedical.com
 info@cobramedical.com
 België / Belgique
 T 070 66 05 70
 F 070 66 05 72

L-GND-SU-01-ENG Ver.4.0 (2015-08)

Nasal Biliary Drainages

Single Use

Nasal Biliary Drainages


- Various Types and Sizes Available
Fit for all anatomic cases
- Well Balanced Wall Diameter
Maximum inner lumen with highest form stability & kink resistance
- Tapered & Round Ends
Atraumatic implantation



Manufactured by:

Medi-Globe GmbH | Medi-Globe-Straße 1-5 | 83101 Achenmühle, Germany
 ☎ +49 8032/973-379 | 📠 +49 8032/973-339 | ✉ sales@medi-globe.de


Nasal Biliary Drainages



Nasal Biliary Drainages Sets

Components: Drainage Catheter with Position Marker, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube Length	Connecting Tube Length
		French Size	Length		
GND-02-05-290S	Standard	5.0	250 cm	50 cm	135 cm
GND-02-06-290S	Standard	6.0	250 cm	50 cm	135 cm
GND-02-07-290S	Standard	7.0	250 cm	50 cm	135 cm
GND-02-10-290S	Standard	10.0	250 cm	50 cm	135 cm



Nasal Biliary Drainages Sets

Components: Drainage Catheter with Position Marker, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube Length	Connecting Tube Length
		French Size	Length		
GND-03-05-290S	Standard with Pigtail	5.0	250 cm	50 cm	135 cm
GND-03-06-290S	Standard with Pigtail	6.0	250 cm	50 cm	135 cm
GND-03-07-290S	Standard with Pigtail	7.0	250 cm	50 cm	135 cm
GND-03-10-290S	Standard with Pigtail	10.0	250 cm	50 cm	135 cm



Nasal Biliary Drainages Sets

Components: Drainage Catheter, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube Length	Connecting Tube Length
		French Size	Length		
GND-01-05-290S	Pigtail	5.0	250 cm	50 cm	135 cm
GND-01-06-290S	Pigtail	6.0	250 cm	50 cm	135 cm
GND-01-07-290S	Pigtail	7.0	250 cm	50 cm	135 cm
GND-01-10-290S	Pigtail	10.0	250 cm	50 cm	135 cm



Nasal Biliary Drainages Sets

Components: Drainage Catheter, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube Length	Connecting Tube Length
		French Size	Length		
GND-04-05-290S	O – Type	5.0	250 cm	50 cm	135 cm
GND-04-06-290S	O – Type	6.0	250 cm	50 cm	135 cm
GND-04-07-290S	O – Type	7.0	250 cm	50 cm	135 cm



Nasal Biliary Drainages Sets

Components: Drainage Catheter, Nasal Tube, Connecting Tube and Multi-Medfix Injection Connector

Catalog No.	Catheter Shape	Drainage Catheter		Nasal Tube Length	Connecting Tube Length
		French Size	Length		
GND-05-05-290S	J – Type	5.0	250 cm	50 cm	135 cm
GND-05-06-290S	J – Type	6.0	250 cm	50 cm	135 cm
GND-05-07-290S	J – Type	7.0	250 cm	50 cm	135 cm

All Drainage Catheters are also available individually (Catalog set no. without “S”).