

MediLime Tummy Button

12 french

3ml ballon

0320-12-0.8 MediLime Tummy Button 12 fr 0,8 cm
 0320-12-1.0 MediLime Tummy Button 12 fr 1,0 cm
 0320-12-1.2 MediLime Tummy Button 12 fr 1,2 cm
 0320-12-1.5 MediLime Tummy Button 12 fr 1,5 cm
 0320-12-1.7 MediLime Tummy Button 12 fr 1,7 cm
 0320-12-2.0 MediLime Tummy Button 12 fr 2,0 cm
 0320-12-2.3 MediLime Tummy Button 12 fr 2,3 cm
 0320-12-2.5 MediLime Tummy Button 12 fr 2,5 cm
 0320-12-2.7 MediLime Tummy Button 12 fr 2,7 cm
 0320-12-3.0 MediLime Tummy Button 12 fr 3,0 cm
 0320-12-3.5 MediLime Tummy Button 12 fr 3,5 cm
 0320-12-4.0 MediLime Tummy Button 12 fr 4,0 cm

14 french

4ml ballon

0320-14-0.8 MediLime Tummy Button 14 fr 0,8 cm
 0320-14-1.0 MediLime Tummy Button 14 fr 1,0 cm
 0320-14-1.2 MediLime Tummy Button 14 fr 1,2 cm
 0320-14-1.5 MediLime Tummy Button 14 fr 1,5 cm
 0320-14-1.7 MediLime Tummy Button 14 fr 1,7 cm
 0320-14-2.0 MediLime Tummy Button 14 fr 2,0 cm
 0320-14-2.3 MediLime Tummy Button 14 fr 2,3 cm
 0320-14-2.5 MediLime Tummy Button 14 fr 2,5 cm
 0320-14-2.7 MediLime Tummy Button 14 fr 2,7 cm
 0320-14-3.0 MediLime Tummy Button 14 fr 3,0 cm
 0320-14-3.5 MediLime Tummy Button 14 fr 3,5 cm
 0320-14-4.0 MediLime Tummy Button 14 fr 4,0 cm
 0320-14-4.5 MediLime Tummy Button 14 fr 4,5 cm
 0320-14-5.0 MediLime Tummy Button 14 fr 5,0 cm

16 french

5ml ballon

0320-16-0.8 MediLime Tummy Button 16 fr 0,8 cm
 0320-16-1.0 MediLime Tummy Button 16 fr 1,0 cm
 0320-16-1.2 MediLime Tummy Button 16 fr 1,2 cm
 0320-16-1.5 MediLime Tummy Button 16 fr 1,5 cm
 0320-16-1.7 MediLime Tummy Button 16 fr 1,7 cm
 0320-16-2.0 MediLime Tummy Button 16 fr 2,0 cm
 0320-16-2.3 MediLime Tummy Button 16 fr 2,3 cm
 0320-16-2.5 MediLime Tummy Button 16 fr 2,5 cm
 0320-16-2.7 MediLime Tummy Button 16 fr 2,7 cm
 0320-16-3.0 MediLime Tummy Button 16 fr 3,0 cm
 0320-16-3.5 MediLime Tummy Button 16 fr 3,5 cm
 0320-16-4.0 MediLime Tummy Button 16 fr 4,0 cm
 0320-16-4.5 MediLime Tummy Button 16 fr 4,5 cm
 0320-16-5.0 MediLime Tummy Button 16 fr 5,0 cm

Inhoud verpakking MediLime Tummy Button:

| | |
|------------------------|-----------|
| MediLime Tummy Button | Spuit 6cc |
| Invoergeleider | Gaasje |
| Lubricating gel 5 gram | |

Productkenmerken

| | |
|---------------------|-------------------------------|
| Materiaal | Siliconen |
| DEHP-vrij | JA |
| Latex vrij | JA |
| Disposabel | JA |
| Steriel | JA |
| Sterilisatiemethode | EO |
| Zichtbaarheid | Rx |
| Verpakking | Omdoos met blister verpakking |
| Type Barcode | GSI |
| Eenheid | Per stuk |
| CE-markering | CE0473 |

18 french

5ml ballon

0320-18-1.0 MediLime Tummy Button 18 fr 1,0 cm
 0320-18-1.2 MediLime Tummy Button 18 fr 1,2 cm
 0320-18-1.5 MediLime Tummy Button 18 fr 1,5 cm
 0320-18-1.7 MediLime Tummy Button 18 fr 1,7 cm
 0320-18-2.0 MediLime Tummy Button 18 fr 2,0 cm
 0320-18-2.3 MediLime Tummy Button 18 fr 2,3 cm
 0320-18-2.5 MediLime Tummy Button 18 fr 2,5 cm
 0320-18-2.7 MediLime Tummy Button 18 fr 2,7 cm
 0320-18-3.0 MediLime Tummy Button 18 fr 3,0 cm
 0320-18-3.5 MediLime Tummy Button 18 fr 3,5 cm
 0320-18-4.0 MediLime Tummy Button 18 fr 4,0 cm
 0320-18-4.5 MediLime Tummy Button 18 fr 4,5 cm
 0320-18-5.0 MediLime Tummy Button 18 fr 5,0 cm

20 french

8ml ballon

0320-20-0.8 MediLime Tummy Button 20 fr 0,8 cm
 0320-20-1.0 MediLime Tummy Button 20 fr 1,0 cm
 0320-20-1.2 MediLime Tummy Button 20 fr 1,2 cm
 0320-20-1.5 MediLime Tummy Button 20 fr 1,5 cm
 0320-20-1.7 MediLime Tummy Button 20 fr 1,7 cm
 0320-20-2.0 MediLime Tummy Button 20 fr 2,0 cm
 0320-20-2.3 MediLime Tummy Button 20 fr 2,3 cm
 0320-20-2.5 MediLime Tummy Button 20 fr 2,5 cm
 0320-20-2.7 MediLime Tummy Button 20 fr 2,7 cm
 0320-20-3.0 MediLime Tummy Button 20 fr 3,0 cm
 0320-20-3.5 MediLime Tummy Button 20 fr 3,5 cm
 0320-20-4.0 MediLime Tummy Button 20 fr 4,0 cm
 0320-20-4.5 MediLime Tummy Button 20 fr 4,5 cm
 0320-20-5.0 MediLime Tummy Button 20 fr 5,0 cm

24 french

8ml ballon

0320-24-1.5 MediLime Tummy Button 24 fr 1,5 cm
 0320-24-1.7 MediLime Tummy Button 24 fr 1,7 cm
 0320-24-2.0 MediLime Tummy Button 24 fr 2,0 cm
 0320-24-2.3 MediLime Tummy Button 24 fr 2,3 cm
 0320-24-2.5 MediLime Tummy Button 24 fr 2,5 cm
 0320-24-2.7 MediLime Tummy Button 24 fr 2,7 cm
 0320-24-3.0 MediLime Tummy Button 24 fr 3,0 cm
 0320-24-3.5 MediLime Tummy Button 24 fr 3,5 cm
 0320-24-4.0 MediLime Tummy Button 24 fr 4,0 cm
 0320-24-4.5 MediLime Tummy Button 24 fr 4,5 cm
 0320-24-5.0 MediLime Tummy Button 24 fr 5,0 cm

