

Stone Extraction Balloons

2 & 3 Lumen – Single Use

Stone Extraction Balloons

➤ Version in  Latex-Free Balloon Quality
Eliminates patient risk for allergic reactions

➤ Multi-Stage Extractor
*Flexible diameter adjustment from
9 - 16 mm in One Balloon*

*Including a set of 3 syringes for predefined
inflation diameters of 9 mm, 13 mm & 16 mm*

➤ Versions with Contrast Injection **above** or
below the Balloon



 **Medi-Globe**[®]
The Spirit of Care.[™]

Stone Extraction Balloon

- Version in Latex-Free Balloon Quality
Eliminates patient risk for allergic reactions
- Radiopaque Marking at the Distal and Proximal End of the Balloon
Precise positioning
- One Centimeter Marking on the Catheter
Macroscopically controlled penetration depth

Multi-Stage Extractor

additionally:

- Flexible Diameter Adjustment from 9 - 16 mm in One Balloon
Time-saving and cost-effective
- Including Set of 3 Syringes for Predefined Inflation of 9 mm, 13 mm & 16 mm
Safe for exact inflation



Contrast Injection Above the Balloon

3-Lumen Stone Extraction Balloon

Catalog No.	Balloon Material	Balloon Diameter	Air Inflation	Liquid Inflation	Catheter Diameter	Working Length	Guide Wire Compatibility
GEB-02-25-035		16 mm	5.0 ml	2.5 ml	7.5 Fr.	200 cm	.035"
GEB-06-25-035		16 mm	5.0 ml	2.5 ml	7.5 Fr.	200 cm	.035"

3-Lumen Multi-Stage Extractor, including a Set of 3 Syringes for Predefined Inflation Diameter of 9 mm, 13 mm & 16 mm

Catalog No.	Balloon Material	Balloon Diameter	Air Inflation	Liquid Inflation	Catheter Diameter	Working Length	Guide Wire Compatibility
GEB-03-25-035		9 - 13 - 16 mm	1.5 ml – 9 mm 2.5 ml – 13 mm 5.0 ml – 16 mm	Non-applicable	7.5 Fr.	200 cm	.035"
GEB-07-25-035		9 - 13 - 16 mm	1.5 ml – 9 mm 2.5 ml – 13 mm 5.0 ml – 16 mm	Non-applicable	7.5 Fr.	200 cm	.035"



2-Lumen Stone Extraction Balloon

Catalog No.	Balloon Material	Balloon Diameter	Air Inflation	Liquid Inflation	Catheter Diameter	Working Length	Guide Wire Compatibility
GEB-01-17-025		12 mm	3.0 ml	1.5 ml	5 Fr.	200 cm	.025"
GEB-05-17-025		12 mm	3.0 ml	1.5 ml	5 Fr.	200 cm	.025"
GEB-01-23-035		16 mm	5.0 ml	2.5 ml	7 Fr.	200 cm	.035"
GEB-05-23-035		16 mm	5.0 ml	2.5 ml	7 Fr.	200 cm	.035"





Contrast Injection below the Balloon

3-Lumen Stone Extraction Balloon

Catalog No.	Balloon Material	Balloon Diameter	Air Inflation	Liquid Inflation	Catheter Diameter	Working Length	Guide Wire Compatibility
GEB-02-25-035P		16 mm	5,0 ml	2,5 ml	7,5 Fr.	200 cm	.035"
GEB-06-25-035P		16 mm	5,0 ml	2,5 ml	7,5 Fr.	200 cm	.035"

3-Lumen Multi-Stage Extractor, including a Set of 3 Syringes for Predefined Inflation Diameter of 9 mm, 13 mm & 16 mm


Catalog No.	Balloon Material	Balloon Diameter	Air Inflation	Liquid Inflation	Catheter Diameter	Working Length	Guide Wire Compatibility
GEB-03-25-035P		9 - 13 - 16 mm	1,5 ml – 9 mm 2,5 ml – 13 mm 5,0 ml – 16 mm	Non-applicable	7.5 Fr.	200 cm	.035"
GEB-07-25-035P		9 - 13 - 16 mm	1,5 ml – 9 mm 2,5 ml – 13 mm 5,0 ml – 16 mm	Non-applicable	7.5 Fr.	200 cm	.035"

Note for Stone Extraction Balloon Versions with Latex:

Caution: This Product Contains Latex Which May Cause Allergic Reactions.

Stone Extraction Balloons

2 & 3 Lumen – Single Use



Ver. 7.0 L-GEB-SU-01-ENG (2015-07)