


Radifocus® Introducer II Standard kits (1/2)

Introducers



Kit of vascular introducer with haemostatic valve with excellent blood tightness, easy and non traumatic vessel insertion, patented valve design for a unique features package.

- **Introducer with extra thin radiopaque wall**
- **Color coded and French size indication**
- **Non traumatic dilator with perfect adjustment to the sheath and safety snap-on connector system**
- **Silicone single wall haemostatic valve with Terumo patented “cross-cut” design**
- **Terumo Nitinol super elastic mini Guide wire available with double end tip (straight or J curved) in A and B kits or spring metallic wire in B and R kits**

Specifications

Kit type	Code number	Inner diameter	Sheath Length	Compatible Guide wire	Guide wire type	Entry venous Catheter	Syringe	Color code
A kit in tray with plastic Guide wire	RS*A40G10SQ	4 Fr	10 cm	0.025" / 0.64 mm	45 cm Straight	20 G x 2" / 0.9 x 51 mm	Yes	Red
	RS*A40K10AQ	4 Fr	10 cm	0.035" / 0.89 mm	45 cm Angled	18 G x 2 1/2" / 1.2 x 64 mm	Yes	Red
	RS*A50K10SQ	5 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight	18 G x 2 1/2" / 1.2 x 64 mm	Yes	Grey
	RS*A60K10SQ	6 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight	18 G x 2 1/2" / 1.2 x 64 mm	Yes	Light green
	RS*A70K10SQ	7 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight	18 G x 2 1/2" / 1.2 x 64 mm	Yes	Orange
	RS*A80K10SQ	8 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight	18 G x 2 1/2" / 1.2 x 64 mm	Yes	Light blue
	RS*A90K10SQ	9 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight	18 G x 2 1/2" / 1.2 x 64 mm	Yes	Black
	RS*A10K10SQ	10 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight	18 G x 2 1/2" / 1.2 x 64 mm	Yes	Purple
RS*A11K10SQ	11 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight	18 G x 2 1/2" / 1.2 x 64 mm	Yes	Yellow	
B kit in tray with plastic Guide wire	RS*B40G10SQ	4 Fr	10 cm	0.025" / 0.64 mm	45 cm Straight	No	No	Red
	RS*B40K10AQ	4 Fr	10 cm	0.035" / 0.89 mm	45 cm Angled	No	No	Red
	RS*B40K10SQ	4 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight	No	No	Red
	RS*B50N10AQ	5 Fr	10 cm	0.038" / 0.97 mm	45 cm Angled	No	No	Grey
	RS*B50N10SQ	5 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight	No	No	Grey
	RS*B50N25AQ	5 Fr	25 cm	0.038" / 0.97 mm	80 cm Angled	No	No	Grey
	RS*B60N10SQ	6 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight	No	No	Light green
	RS*B60N25AQ	6 Fr	25 cm	0.038" / 0.97 mm	80 cm Angled	No	No	Light green
	RS*B70N10SQ	7 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight	No	No	Orange
	RS*B70N25AQ	7 Fr	25 cm	0.038" / 0.97 mm	80 cm Angled	No	No	Orange
	RS*B80N10SQ	8 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight	No	No	Light blue
	RS*B80N25AQ	8 Fr	25 cm	0.038" / 0.97 mm	80 cm Angled	No	No	Light blue
	RS*B90N10SQ	9 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight	No	No	Black
	RS*B90N25AQ	9 Fr	25 cm	0.038" / 0.97 mm	80 cm Angled	No	No	Black
	RS*B10N10SQ	10 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight	No	No	Purple
	RS*B10N25AQ	10 Fr	25 cm	0.038" / 0.97 mm	80 cm Angled	No	No	Purple
RS*B11N10SQ	11 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight	No	No	Yellow	
RS*B11N25AQ	11 Fr	25 cm	0.038" / 0.97 mm	80 cm Angled	No	No	Yellow	

A Kit includes: introducer sheath, dilator, mini Guide wire, entry venous catheter, 2.5 cc syringe

B Kit includes: introducer sheath, dilator, mini Guide wire

C Kit includes: introducer sheath, dilator

R Kit includes: introducer sheath, dilator, mini Guide wire, 18 G x 2 3/4" / 1.2 x 70 mm entry metallic needle

Other code numbers are available on special demand. For any further information, please contact your Terumo local representative

Units per box: 5 for kits in tray and 10 for kits in pouch

Radifocus® Introducer II Standard kits (2/2)

Introducers

Specifications

Kit type	Code number	Inner diameter	Sheath Length	Compatible Guide wire	Guide wire type	Entry venous Catheter	Syringe	Color code
B kit in tray with spring double end Guide wire	RS*B40K10MQ	4 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Red
	RS*B50K10MQ	5 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Grey
	RS*B50N10MQ	5 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Grey
	RS*B60K10MQ	6 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Light green
	RS*B60N10MQ	6 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Light green
	RS*B70K10MQ	7 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Orange
	RS*B70N10MQ	7 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Orange
	RS*B80K10MQ	8 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Light blue
	RS*B80N10MQ	8 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Light blue
	RS*B90K10MQ	9 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Black
	RS*B90N10MQ	9 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Black
	RS*B10K10MQ	10 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Purple
	RS*B10N10MQ	10 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Purple
RS*B11K10MQ	11 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Yellow	
RS*B11N10MQ	11 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Yellow	
B kit in pouch with spring double end Guide wire	RS*B40K10MR	4 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	No	No	Red
	RS*B50N10MRD	5 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Grey
	RS*B60N10MRD	6 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Light green
	RS*B70N10MRD	7 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Orange
	RS*B80N10MRD	8 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Light blue
	RS*B90N10MRD	9 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Black
	RS*B10N10MRD	10 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Purple
	RS*B11N10MRD	11 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	No	No	Yellow
C Kit in pouch without Guide wire	RS*C40G10NR	4 Fr	10 cm	0.025" / 0.64 mm	No mini Guide wire	No	No	Red
	RS*C40K10NR	4 Fr	10 cm	0.035" / 0.89 mm	No mini Guide wire	No	No	Red
	RS*C40K25NR	4 Fr	25 cm	0.035" / 0.89 mm	No mini Guide wire	No	No	Red
	RS*C50N10NR	5 Fr	10 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Grey
	RS*C60N10NR	6 Fr	10 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Light green
	RS*C60N25NR	6 Fr	25 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Light green
	RS*C70N10NR	7 Fr	10 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Orange
	RS*C70N25NR	7 Fr	25 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Orange
	RS*C80N10NR	8 Fr	10 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Light blue
	RS*C80N25NR	8 Fr	25 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Light blue
	RS*C90N10NR	9 Fr	10 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Black
	RS*C10N10NR	10 Fr	10 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Purple
RS*C11N10NR	11 Fr	10 cm	0.038" / 0.97 mm	No mini Guide wire	No	No	Yellow	
R kit in tray with spring double end Guide wire + metallic needle	RS*R40K10MQ	4 Fr	10 cm	0.035" / 0.89 mm	45 cm Straight & J	18 G x 2 3/4" / 1.2 x 70 mm	No	Red
	RS*R50N10MQ	5 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	18 G x 2 3/4" / 1.2 x 70 mm	No	Grey
	RS*R60N10MQ	6 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	18 G x 2 3/4" / 1.2 x 70 mm	No	Light green
	RS*R70N10MQ	7 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	18 G x 2 3/4" / 1.2 x 70 mm	No	Orange
	RS*R80N10MQ	8 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	18 G x 2 3/4" / 1.2 x 70 mm	No	Light blue
	RS*R90N10MQ	9 Fr	10 cm	0.038" / 0.97 mm	45 cm Straight & J	18 G x 2 3/4" / 1.2 x 70 mm	No	Black

A Kit includes: introducer sheath, dilator, mini Guide wire, entry venous catheter, 2.5 cc syringe

B Kit includes: introducer sheath, dilator, mini Guide wire

C Kit includes: introducer sheath, dilator

R Kit includes: introducer sheath, dilator, mini Guide wire, 18 G x 2 3/4" / 1.2 x 70 mm entry metallic needle

Other code numbers are available on special demand. For any further information, please contact your Terumo local representative

Units per box: 5 for kits in tray and 10 for kits in pouch

Radifocus® Introducer II with Radiopaque marker

Introducers

Kit of vascular introducer with haemostatic valve with excellent blood tightness, easy and non traumatic vessel insertion, patented valve design for a unique features package and with distal radiopaque marker.

- Introducer with distal highly visible Platinum/Iridium band incorporated in extra thin wall.
- Color coded and French size indication
- Non traumatic dilator with perfect adjustment to the sheath and safety snap-on connector system
- Silicone single wall haemostatic valve with Terumo patented “cross-cut” design

Specifications

Kit type	Code number	Inner diam.	Sheath length	Compatible Guide wire	Guide wire type	Distal marker	Color code
C Kit in pouch	RB*C6SN06NC	6 Fr	6 cm	0.038" / 0.97 mm		Pt/Ir	Light green
	RB*C6SN10NC	6 Fr	10 cm	0.038" / 0.97 mm		Pt/Ir	Light green
	RB*C6SN25NC	6 Fr	25 cm	0.038" / 0.97 mm		Pt/Ir	Light green
	RB*C7SN06NC	7 Fr	6 cm	0.038" / 0.97 mm		Pt/Ir	Orange
	RB*C7SN10NC	7 Fr	10 cm	0.038" / 0.97 mm		Pt/Ir	Orange
	RB*C7SN25NC	7 Fr	25 cm	0.038" / 0.97 mm		Pt/Ir	Orange
	RB*C8SN06NC	8 Fr	6 cm	0.038" / 0.97 mm		Pt/Ir	Light blue
	RB*C8SN10NC	8 Fr	10 cm	0.038" / 0.97 mm		Pt/Ir	Light blue
	RB*C8SN25NC	8 Fr	25 cm	0.038" / 0.97 mm		Pt/Ir	Light blue

C Kit includes: introducer sheath, dilator

Units per box: 10

Radifocus® Introducer II Pediatric kits

Introducers

Kit of vascular introducer with haemostatic valve with excellent blood tightness, easy and non traumatic vessel insertion, patented valve design for a unique features package.

- **Introducer with extra thin radiopaque wall**
- **Color coded and French size indication**
- **Non traumatic dilator with perfect adjustment to the sheath and safety snap-on connector system**
- **Silicone single wall haemostatic valve with Terumo patented “cross-cut” design**
- **Terumo Nitinol super elastic mini-guidewire available with double end tips (straight or J curved)**

Specifications

Kit Type	Code number	Inner diameter	Sheath length	Compatible Guide wire	Guide wire type	Entry venous catheter	Syringe	Color Code
A Kit in tray Plastic Guide wire	RS*A40G05SQ	4 Fr	5cm	0.025" / 0.64 mm	45 cm straight	20 G x 2" / 0.9 x 51 mm	Yes	Red
	RS*A40G07SQ	4 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	20 G x 2" / 0.9 x 51 mm	Yes	Red
	RS*A50G05SQ	5 Fr	5 cm	0.025" / 0.64 mm	45 cm straight	20 G x 2" / 0.9 x 51 mm	Yes	Grey
	RS*A50G07SQ	5 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	20 G x 2" / 0.9 x 51 mm	Yes	Grey
	RS*A60G05SQ	6 Fr	5 cm	0.025" / 0.64 mm	45 cm straight	20 G x 2" / 0.9 x 51 mm	Yes	Light green
	RS*A60G07SQ	6 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	20 G x 2" / 0.9 x 51 mm	Yes	Light green
B Kit in tray Plastic Guide wire	RS*B40G07SQ	4 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	No	No	Red
	RS*B50G07SQ	5 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	No	No	Grey
	RS*B60G07SQ	6 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	No	No	Light green

A Kit includes: introducer sheath, dilator, mini Guide wire, entry venous catheter, 2.5cc syringe

B Kit includes: introducer sheath, dilator, mini Guide wire

Units per box: 5

Radifocus® Introducer II Transradial kits

Introducers



Standard short sheath

Introducer kit II concentrates Terumo's minimally invasive technology know-how:

- Specifically for radial access, dilators are tapered on each guidewire to minimise vessel trauma
- Smooth sheath to dilator transition
- Version with metal needle and spring wire

Specifications

Kit in tray with metal needle and spring mini Guide wire

Code number	Inner sheath diameter	Sheath Length	Compatible Guide wire	Mini Guide wire type	Metal needle	Syringe	Color code
RT-R40A07PQ	4 Fr	7 cm	0.018" / 0.46 mm	45 cm straight	22 G x 1 ^{1/2} " / 0,7 x 38 mm	no	Red
RT-R40A10PQ	4 Fr	10 cm	0.018" / 0.46 mm	45 cm straight	22 G x 1 ^{1/2} " / 0,7 x 38 mm	no	Red
RT-R40D10PQ	4 Fr	10 cm	0.021" / 0.53 mm	45 cm straight	21 G x 1 ^{1/2} " / 0,8 x 38 mm	no	Red
RT-R40G07PQ	4 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	20 G x 1 ^{1/2} " / 0,9 x 38 mm	no	Red
RT-R40G10PQ	4 Fr	10 cm	0.025" / 0.64 mm	45 cm straight	20 G x 1 ^{1/2} " / 0,9 x 38 mm	no	Red
RT-R50A07PQ	5 Fr	7 cm	0.018" / 0.46 mm	45 cm straight	22 G x 1 ^{1/2} " / 0,7 x 38 mm	no	Grey
RT-R50A10PQ	5 Fr	10 cm	0.018" / 0.46 mm	45 cm straight	22 G x 1 ^{1/2} " / 0,7 x 38 mm	no	Grey
RT-R50D10PQ	5 Fr	10 cm	0.021" / 0.53 mm	45 cm straight	21 G x 1 ^{1/2} " / 0,8 x 38 mm	no	Grey
RT-R50G07PQ	5 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	20 G x 1 ^{1/2} " / 0,9 x 38 mm	no	Grey
RT-R50G10PQ	5 Fr	10 cm	0.025" / 0.64 mm	45 cm straight	20 G x 1 ^{1/2} " / 0,9 x 38 mm	no	Grey
RT-R60A07PQ	6 Fr	7 cm	0.018" / 0.46 mm	45 cm straight	22 G x 1 ^{1/2} " / 0,7 x 38 mm	no	Light green
RT-R60A10PQ	6 Fr	10 cm	0.018" / 0.46 mm	45 cm straight	22 G x 1 ^{1/2} " / 0,7 x 38 mm	no	Light green
RT-R60D10PQ	6 Fr	10 cm	0.021" / 0.53 mm	45 cm straight	21 G x 1 ^{1/2} " / 0,8 x 38 mm	no	Light green
RT-R60G07PQ	6 Fr	7 cm	0.025" / 0.64 mm	45 cm straight	20 G x 1 ^{1/2} " / 0,9 x 38 mm	no	Light green
RT-R60G10PQ	6 Fr	10 cm	0.025" / 0.64 mm	45 cm straight	20 G x 1 ^{1/2} " / 0,9 x 38 mm	no	Light green
RT-R70D10PQ	7 Fr	10 cm	0.021" / 0.53 mm	45 cm straight	21 G x 1 ^{1/2} " / 0,8 x 38 mm	no	Orange

Units per box: 5

Radifocus® Introducer II M Coat™ Hydrophilic Radial Access Introducer kit

Introducers



M Coat™ long hydrophilic radial sheath

Terumo's M Coat™ sheath offers smoother insertion and removal and minimizes the problem of spasm, making it ideal for routine TRI* procedures. Retrieving the M Coat™ sheath from the radial artery requires only half the force compared to conventional uncoated sheaths**. This helps minimizing patient discomfort.

* Trans Radial Intervention
** Confer to Tractors study

Specifications

Code number	Inner sheath diameter	Sheath length	Compatible Guide wire	Mini Guide wire type	Entry venous catheter	Syringe	Color code
A Kit in tray Plastic Guide wire							
RM*AF5J10SQW	5 Fr	10 cm	0.025" / 0.64 mm	45 cm straight	20 G x 2"/ 0,9 x 51 mm	Yes	Grey
RM*AF5J16SQW	5 Fr	16 cm	0.025" / 0.64 mm	80 cm straight	20 G x 2"/ 0,9 x 51 mm	Yes	Grey
RM*AF5J25SQW	5 Fr	25 cm	0.025" / 0.64 mm	80 cm straight	20 G x 2"/ 0,9 x 51 mm	Yes	Grey
RM*AF6G10HQW	6 Fr	10 cm	0.025" / 0.64 mm	45 cm short angle	20 G x 1 ^{1/4} / 0,9 x 32 mm	Yes	Light green
RM*AF6J10SQW	6 Fr	10 cm	0.025" / 0.64 mm	45 cm straight	20 G x 2"/ 0,9 x 51 mm	Yes	Light green
RM*AF6J16SQW	6 Fr	16 cm	0.025" / 0.64 mm	80 cm straight	20 G x 2"/ 0,9 x 51 mm	Yes	Light green
RM*AF6J25SQW	6 Fr	25 cm	0.025" / 0.64 mm	80 cm straight	20 G x 2"/ 0,9 x 51 mm	Yes	Light green
B Kit in tray spring wire							
RM*BF5F10PQ	5 Fr	10 cm	0.021" / 0.53 mm	45 cm straight	No	No	Grey
RM*BF5J10PQ	5 Fr	10 cm	0.025" / 0.64 mm	45 cm straight	No	No	Grey
RM*BF5F16PQ	5 Fr	16 cm	0.021" / 0.53 mm	80 cm straight	No	No	Grey
RM*BF5J16PQ	5 Fr	16 cm	0.025" / 0.64 mm	80 cm straight	No	No	Grey
RM*BF5F25PQ	5 Fr	25 cm	0.021" / 0.53 mm	80 cm straight	No	No	Grey
RM*BF5J25PQ	5 Fr	25 cm	0.025" / 0.64 mm	80 cm straight	No	No	Grey
RM*BF6F10PQ	6 Fr	10 cm	0.021" / 0.53 mm	45 cm straight	No	No	Light green
RM*BF6J10PQ	6 Fr	10 cm	0.025" / 0.64 mm	45 cm straight	No	No	Light green
RM*BF6F16PQ	6 Fr	16 cm	0.021" / 0.53 mm	80 cm straight	No	No	Light green
RM*BF6J16PQ	6 Fr	16 cm	0.025" / 0.64 mm	80 cm straight	No	No	Light green
RM*BF6F25PQ	6 Fr	25 cm	0.021" / 0.53 mm	80 cm straight	No	No	Light green
RM*BF6J25PQ	6 Fr	25 cm	0.025" / 0.64 mm	80 cm straight	No	No	Light green

Units per box: 5